

## Patient Specific Tubular Slide Sheet



interwe

healthcare



## Care

Our patient specific slide sheets are designed to be kept with the same patient for the duration of their care and discarded once their usage has finished.

Though unsuitable for laundering, it may be possible to use a diluted cleaning solution for lightly soiled sheets.

## Usage

Slide sheets are not intended for lifting patients, all repositioning manoeuvres should be accomplished following safe handling practices.

Please inspect items prior to use and be mindful of normal wear and tear over time when slide sheets should be replaced as deemed necessary.

If in any doubt, consult your Health and Safety or Moving and Handling advisor before use.