


# Elite Reusable Tubular Slide Sheet

**interweave**  
healthcare

<b>Product Code:</b>	PH407-90-90---TEAA77
<b>NHS Supply Chain:</b>	FES10784
<b>Size:</b>	90cm x 90cm
<b>Colour:</b>	Teal
<b>Material:</b>	100% silicone coated polyester
<b>Weights:</b>	127g / 70gsm (+/- 5%)
<b>Description:</b>	ELITE SLIDETEX REUSABLE TUBULAR SLIDE SHEET 90X90#CM TEAL



## Care

Our reusable slide sheets can be laundered following the standard hospital laundry guidelines HSG (95)18 and are machine washable up to 71 degrees. Drip dry, or if necessary, a low temperature tumble dry may be used with care.

Tumble dried slide sheets will likely develop a wrinkled appearance but this should not affect performance.

Alternatively, sheets may be wiped clean using a dilute chlorine-based solution or a good quality sterilisation fluid.

## Usage

Slide sheets are not intended for lifting patients, all repositioning manoeuvres should be accomplished following safe handling practices.

Please inspect items prior to use and be mindful of normal wear and tear over time when slide sheets should be replaced as deemed necessary.

If in doubt, always consult your Health and Safety or Moving and Handling advisor before use.